

P-AMINOBENZOIC ACID

 INGREDIENTS
 CAS NO
 %
 8HR OEL

 p-aminobenzoic acid
 150-13-0
 >=99

UN No: Not Applicable
Hazchem Code: Not Applicable
DG Class: Not Applicable
Subsidiary Risk: Not Applicable
Packing Group: Not Applicable
Poisons Schedule: Not Applicable



NFPA Rating:

- ▶ 0: Minimum
- ▶ 1: Low
- ▶ 2: Moderate
- ▶ 3: High
- ▶ 4: Extreme

HEALTH HAZARD INFORMATION













Harmful if swallowed. Irritating to eyes.

Acute Health Irritating to respiratory system.

Effects: Irritating to skin.

May cause SENSITISATION by inhalation.
May cause SENSITISATION by skin contact.
Possible risk of irreversible effects.

PRECAUTIONS FOR USE













Appropriate engineering controls:

General Exhaust Ventilation adequate.

Glasses:

Consider chemical goggles.

Respirator:

Particulate. (AS/NZS 1716 & 1715, EN 143:2000 & 149:001, ANSI Z88 or national equivalent)

Storage and Transportation: Store in cool, dry, protected area. Dispose of this material and its container at hazardous or special waste collection point. Keep out of reach of children. Keep away from food, drink and animal feeding stuffs.

Fire/Explosion Hazard:

Toxic smoke/fumes in a fire. Dispose of this material and its container at hazardous or special waste collection point.

PROPERTIES





Solid. Mixes with water.Combustible.

EMERGENCY









FIRST AID

Swallowed: If more than 15 mins from Doctor, INDUCE VOMITING (if conscious).

Eye: Wash with running water.

Skin: Remove contaminated clothing. Wash with soap &

water.

Inhaled: Fresh air. Rest, keep warm. If breath shallow, give

oxygen. Medical attention.

Advice To
Doctor: Assist ventilation Hypotension - Trendelenburg
Position. Methylene blue for severe poisoning.

Fire Fighting: Water spray/ fog.

Eliminate ignition sources. Avoid dust. Sweep shovel to safe place. Dispose of this material and its container at hazardous or special waste collection point. This material and its container must be

SAFE STORAGE WITH OTHER CLASSIFIED CHEMICALS





Explosive Radioactive

Spills and

Disposal:

A

— Must not be stored together